

INSTALLATION, OPERATION AND MAINTENANCE MANUAL

CLIMATE CONTROL DEVICE 4000-01



TYPE CE, CONTINUOUS FLOW, SUPPLIED-AIR RESPIRATOR
THESE RESPIRATORS ARE APPROVED ONLY IN THE FOLLOWING CONFIGURATIONS

RESPIRATOR COMPONENTS																												
TC No.	Protection 1	M O D E L	Alternate Helmets		Alternate Flow Regulator and Belt			Alternate Capes			Alternate Quick Disconnect			Alternate Hoses			Breath- ing Air Tube	Alternate Visor		Alternate Lenses			Cautions and Limitations 2					
			NV2001	02-803	NV2016	4000-01	NV2015	4000-20	NV2002	NV2002L	NV2002XL	02-813	NV2025	NV2024	NV2032	NV2034		NV2035	NV2028	NV2029	NV2035	NV2036		NV2021	NV2004	02-804	NV2018	NV2031
19C-363	SA/CF/SB	NOVA 2000	X		X	X		X	X	X		X	X	X		X	X			X	X		X	X	X			ABCDEJMNOS
19C-417	SA/CF/SB	ASTRO		X	X						X	X	X	X		X	X			X		X				X	X	ABCDEJMNOS
19C-422	SA/CF/SB	NOVA 2000	X				X		X	X	X				X	X		X	X	X	X		X	X	X			ABCDEJMNOS
19C-423	SA/CF/SB	ASTRO		X			X							X	X		X	X	X	X		X				X	X	ABCDEJMNOS

1. PROTECTION

CF- Continuous flow SA- Supplied Air SB- Supplied - Air Abrasive Blast

2. CAUTIONS AND LIMITATIONS

- A-** Not for use in atmospheres containing less than 19.5 percent oxygen.
- B-** Not for use in atmospheres immediately dangerous to life and health.
- C-** Do not exceed maximum use concentrations established by regulatory standards.
- D-** Air-line respirators can be used only when the respirators are supplied with respirable air meeting the requirements of CGA G-7.1 Grade D or Higher Quality.
- E-** Use only the pressure ranges and hose lengths specified in the User's instructions.
- J-** Failure to properly use and maintain this products could result in injury or death.
- M-** All approved respirators shall be selected, fitted, used and maintained in accordance with MSHA, OSHA, and other applicable regulations.
- N-** Never substitute, modify, add, or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer.
- O-** Refer to User's instructions, and/or maintenance manuals for information on use and maintenance of these respirators.
- S-** Special or critical User's instructions and/or specific Limitations apply. Refer to User's instructions.

INTRODUCTION

The 4000-01 Cold Air Tube is NIOSH approved to supply a continuous flow of cool compressed air to the Supplied Air Respirator. All Genuine Parts must be present and correctly assembled to constitute a NIOSH approved respirator.

AIR QUALITY

The respirator must be supplied at all times with clean breathable air Grade D or better. Always use an inlet filter on your compressor - locate your compressor in a clean environment. The AIR QUALITY must be constantly monitored to ensure Grade D air is supplied at all times. A carbon monoxide alarm must always be used.

WARNING

Do not connect the air supply hose to nitrogen, toxic gases, inert gases or any other non Grade D air source.

AIR PRESSURE

The air pressure supplied to the respirator must be continuously monitored to ensure it keeps within the specified range as detailed on the breathing air pressure tables.

WARNING

Make sure the correct table is used. When using 4000-01 use (Fig 3) Breathing Air Pressure Table.

WARNINGS

- **DO NOT** use the 4000-01 until you have been trained in its use, maintenance and limitations by a qualified individual (appointed by your employer) who has extensive knowledge of the 4000-01
- **DO NOT** ingest the ice that may form around the cold air outlet as it could contain oils or toxic element that could cause illness or death.
- **USE ONLY** NIOSH approved NOVA 2000 air supply hose.
- **DO NOT** modify or alter parts of the respirator. The use of non approved parts voids the entire approval of the Respirator.

1. To supply cold air to the respirator connect the outlet to the NV2021 breathing tube and tighten by hand. (Fig 1).
2. Connect the NIOSH approved air supply hose to a breathing air source, supplying Grade D or better quality air. Connect the air supply hose to the 4000-01 quick disconnect tail. (Fig 2).
3. With air flowing through the respirator, adjust the air pressure at the point of attachment. (Fig 2) to the recommended air pressure specified in the BREATHING AIR PRESSURE TABLE (Fig 3).

Fig 1





Fig 2



⚠ WARNING ⚠

Make sure you understand the breathing air pressure table before using the respirator. Make sure you set the correct air pressure for the correct FIG 2 supply hose length.

Fig 3

BREATHING AIR PRESSURE TABLE

AIR SOURCE	BREATHING TUBE ASSEMBLY	AIR SUPPLY HOSE	AIR SUPPLY HOSE LENGTH	MAX NUMBER OF SECTIONS	PRESSURE RANGE (PSIG Air)
Portable or Stationary Compressor	4000-01 Cold Air Tube	NV2028/ NV2029	25	1	55 - 56
			50	1	55 - 57
			100	2	60 - 62
			150	3	65 - 67
			200	4	70 - 72
			250	5	77 - 78
			300	6	80 - 82

4. To obtain cooler air turn the regulator control knob anti clockwise, this will increase the airflow out of the exhaust port (Fig 4). Turn the knob clockwise to obtain air closer to ambient temperature.

⚠ WARNING ⚠


DO NOT USE THE 4000-01 when ambient temperature is below 68; (20;C) as ice could form in the cold air outlet and reduce airflow.

5. With the compressed air flowing into the respirator at the correct pressure range for the correct supply hose length you are now ready to enter the work area. Do not remove your respirator until you are outside the blasting area.


⚠ WARNING ⚠

Leave work area immediately if you are unwell, airflow stops or slows down, breathing becomes difficult, you become dizzy, too cold or if any respirator part is damaged.

Fig 4



CCT 4000 PARTS LIST



Item	PART No.
Top Threaded Cap	4000-03
Air Spinner	4000-04
Body	4000-05
1/4" Quick Disconnect Tail	4000-06
Breaker	4000-07
Regulator Body	4000-08
Spindle & Knob	4000-09
Locking Screw	4000-10
Regulator Assembly	4000-11
Heat Shield	4000-12
Belt	NV2022

⚠ WARNING ⚠

Due to the high noise levels created during abrasive blasting earplugs must be worn at all time. During Abrasive Blasting operations do not use Silica Sand, the use of Silica Sand can cause serious respiratory disease or death, refer: www.osha.gov/Silica.IT69D_1.html



NOTES



NOTES



OTHER PRODUCTS



ISO9001
ACCREDITED COMPANY

DID YOU KNOW

We can also supply you with the following...

NOVA 2000 BLASTING HELMET

The NOVA 2000 has been designed for blasting professionals.

The NOVA 2000 is the most advanced abrasive blasting helmet available, with fully padded pillow foam liner for extra comfort and noise attenuation, wide angle vision, even air distribution assures a pleasant working environment, optional climate control devices. Things cannot get much better.

ASTRO BLASTING HELMET

ASTRO is your number one choice economic Abrasive Blasting Helmet which meets safety standards worldwide.

This is the lightweight, robust helmet you have been waiting for. It helps to reduce fatigue, reduce down time and increase productivity.

HEAD SOCKS & SPRAY HOODS

Stay cleaner with Head Socks and Spray Hoods. One size fits all. The cool comfortable Spray Hood protects from over spray and sanding dust. The Head Sock can be worn as a hygienic liner.

HEAD SOCKS & SPRAY HOODS

The RADEX Airline Filter offers you more capacity, versatility & filtration. You can combine the versatility of either floor or wall mounting with increased filtration capacity.

Our range of optional equipment means you can customise your filter to meet your requirements.



BLASTONE

BLASTONE INTERNATIONAL (AUST) PTY LTD

57 Alexandra Place Murarrie
Queensland 4172 Australia

PO Box 3385 Tingalpa DC
Queensland 4172 Australia

Toll Free 1800 190 190
Tel (07) 3329 5000
Fax (07) 3329 5066
Email sales.au@blastone.com
www.BlastOne.com

AUSTRALIAN OFFICES

Adelaide | Brisbane | Darwin | Mackay |
Melbourne | Perth | Sydney

BLASTONE INTERNATIONAL

4510 Bridgeway Avenue
Columbus Ohio 43219 USA

Toll Free 800 999 1881
Tel (614) 476 3000
Fax (614) 476 3002
Email sales@blastone.com
www.BlastOne.com

NORTH AMERICAN OFFICES

Chicago | Columbus | Los Angeles |
Minneapolis